

Product Datasheet

His Tag antibody (orb1140376)

| | |
|----------------------------|--|
| Description | Rabbit polyclonal antibody to His Tag |
| Species/Host | Rabbit |
| Reactivity | Other |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, ICC, IHC, IP, WB |
| Immunogen | The antiserum was produced by immunizing rabbit with synthesized peptide containing the His epitope (H-H-H-H-H). |
| Form/Appearance | Liquid |
| Concentration | 100ug/100ul |
| Storage | This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt. |
| Note | For research use only |
| Application notes | ELISA 1: 2500; WB 1:500-1000; IP 5 µg/sample; IF 5µg/mL; IHC 5 µg/mL |
| Isotype | IgG |
| Clonality | Polyclonal |
| Dilution Range | ELISA 1: 2500; WB 1:500-1000; IP 5 µg/sample; IF 5µg/mL; IHC 5 µg/mL |
| Expiration Date | 12 months from date of receipt. |