

Product Datasheet

GFP Tag antibody (orb1140378)

Description	Gallus polyclonal antibody to GFP Tag
Species/Host	Gallus
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	The anti-GFP immunogen is a Green Fluorescent Protein (GFP) fusion protein corresponding to the full length amino acid sequence (246aa) derived from the jellyfish <i>Aequorea victoria</i> .
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.