

Product Datasheet

GFP Tag antibody (orb1140400)

Description	Goat polyclonal antibody to GFP Tag
Species/Host	Goat
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC, IP, WB
Immunogen	The anti-GFP immunogen is a Green Fluorescent Protein (GFP) fusion protein corresponding to the full length amino acid sequence (246aa) derived from the jellyfish <i>Aequorea victoria</i> .
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	ELISA 1:40,000; WESTERN BLOT 1:400 - 1:2,000; IMMUNOHISTOCHEMISTRY 1:200 - 1:1,000; IF MICROSCOPY 1:1,000; IP 1:200
Isotype	IgG
Clonality	Polyclonal
Dilution Range	ELISA 1:40,000; WESTERN BLOT 1:400 - 1:2,000; IMMUNOHISTOCHEMISTRY 1:200 - 1:1,000; IF MICROSCOPY 1:1,000; IP 1:200
Expiration Date	12 months from date of receipt.