
Product Datasheet

RFP Tag antibody (orb1140403)

Description

Mouse monoclonal antibody to RFP Tag

Species/Host

Mouse

Reactivity

Other

Conjugation

Unconjugated

Tested Applications

ELISA, IP, WB

Immunogen

RFP from the Discosoma sp. (sea anemone) N-terminal peptide-KLH conjugates.

Form/Appearance

Liquid

Concentration

100ug/100ul

Storage

This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Note

For research use only

Application notes

WB (with ECL): 1:1000-3000 dilution (incubate for one hour at room temperature). IS: 1:500-2000 dilution. For best results with other assays (e.g.: Dot, ELISA, IP, etc), please determine optimal working dilution by titration test.

Isotype

IgG1

Clonality

Monoclonal

Dilution Range

WB (with ECL): 1:1000-3000 dilution (incubate for one hour at room temperature). IS: 1:500-2000 dilution. For best results with other assays (e.g.: Dot, ELISA, IP, etc), please determine optimal working dilution by titration test.

Expiration Date

12 months from date of receipt.