
Product Datasheet

MBP Tag antibody (orb1140414)

Description

Rabbit polyclonal antibody to MBP Tag

Species/Host

Rabbit

Reactivity

Other

Conjugation

Unconjugated

Tested Applications

IHC, IP, WB

Immunogen

MBP epitope tag recombinant protein.

Form/Appearance

Liquid

Concentration

100ug/100ul

Storage

This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Note

For research use only

Application notes

Optimal working dilutions should be determined experimentally by investigators. Prepare working dilution immediately before use. Suggested starting dilutions are as follows. Western Blot: 1:1,000
Immunoprecipitation: 5 µg/mg lysate
Histo/Cytochemistry: 1:200

Isotype

IgG

Clonality

Polyclonal

Dilution Range

Optimal working dilutions should be determined experimentally by investigators. Prepare working dilution immediately before use. Suggested starting dilutions are as follows. Western Blot: 1:1,000
Immunoprecipitation: 5 µg/mg lysate
Histo/Cytochemistry: 1:200

Expiration Date

12 months from date of receipt.