

Product Datasheet

Myc Tag antibody (orb1140416)

Description	Rabbit polyclonal antibody to Myc Tag
Species/Host	Rabbit
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Myc epitope tag peptide E-Q-K-L-I-S-E-E-D-L
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	ELISA 1:1,000 ; WESTERN BLOT 1:500 - 1:1,000 ; IMMUNOHISTOCHEMISTRY User Optimized ; OTHER APPLICATIONS User Optimized
Isotype	IgG
Clonality	Polyclonal
Dilution Range	ELISA 1:1,000 ; WESTERN BLOT 1:500 - 1:1,000 ; IMMUNOHISTOCHEMISTRY User Optimized ; OTHER APPLICATIONS User Optimized
Expiration Date	12 months from date of receipt.