

Product Datasheet

HA Tag antibody (orb1140430)

Description	Mouse monoclonal antibody to HA Tag
Species/Host	Mouse
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IP, WB
Immunogen	A synthetic peptide from influenza hemagglutinin epitope (YPYDVPDYA) coupled to KLH.
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	ICC/IHC/IF/IP 1:200 dilution. WB 1:500-1:1000 dilution.
Isotype	IgG3
Clonality	Monoclonal
Dilution Range	ICC/IHC/IF/IP 1:200 dilution. WB 1:1000 dilution.
Expiration Date	12 months from date of receipt.