

Product Datasheet

His Tag antibody (orb1140432)

Description	Mouse monoclonal antibody to His Tag
Species/Host	Mouse
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IP, WB
Immunogen	6x His synthetic peptide.
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	WB (with ECL): 1:1,000-3,000 dilution (incubate for one hour at room temperature). IS: 1:500-2000 dilution. Other applications are optimized by end users.
Isotype	IgG2b
Clonality	Monoclonal
Dilution Range	WB (with ECL): 1:1,000-3,000 dilution (incubate for one hour at room temperature). IS: 1:500-2000 dilution. Other applications are optimized by end users.
Expiration Date	12 months from date of receipt.