

Product Datasheet

Myc Tag antibody (orb1140433)

Description	Mouse monoclonal antibody to Myc Tag
Species/Host	Mouse
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IP, WB
Immunogen	A synthetic peptide corresponding to residues 410-419 (EQKLISEEDL) coupled to KLH.
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	Western Blot (1:1000), Immunoprecipitation and Immunostaining (1:200). Other dilutions will need optimization by end user.
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	Western Blot (1:1000), Immunoprecipitation and Immunostaining (1:200). Other dilutions will need optimization by end user.
Expiration Date	12 months from date of receipt.