

Product Datasheet

GST Tag antibody (orb1140435)

Description	Mouse monoclonal antibody to GST Tag
Species/Host	Mouse
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	Anti-GST antibody has been raised against whole GST protein.
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is a minimum of one (1) year from date of receipt when stored appropriately, please refer to CoA for lot specific expiration dates.
Note	For research use only
Application notes	Western blot (1:1000), Immunoprecipitation and Immunostaining (1:200).
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	Western blot (1:1000), Immunoprecipitation and Immunostaining (1:200).
Expiration Date	12 months from date of receipt.