

## Product Datasheet

### ATI 2341 (orb1140485)

<b>Description</b>	CXCR4 allosteric agonist; Peptides.
<b>Form/Appearance</b>	Freeze dried solid
<b>Storage</b>	Store dry, frozen and in the dark
<b>Note</b>	For research use only
<b>Protein Sequence</b>	Pal-MGYQKKLRSMTDKYRL
<b>Formula</b>	$C_{104}H_{178}N_{26}O_{25}S_2$
<b>Purity</b>	> 95% by hplc
<b>MW</b>	2256.82 Da
<b>CAS Number</b>	1337878-62-2
<b>Solubility (25°C)</b>	Soluble to 1 mg/ml in water
<b>Expiration Date</b>	12 months from date of receipt.