

## Product Datasheet

### Hexa-D-arginine (orb1140541)

<b>Description</b>	Furin inhibitor; Peptides.
<b>Form/Appearance</b>	Freeze dried solid
<b>Storage</b>	Store dry, frozen and in the dark
<b>Note</b>	For research use only
<b>Protein Sequence</b>	rrrrrr-NH <sub>2</sub>
<b>Formula</b>	C <sub>36</sub> H <sub>75</sub> N <sub>25</sub> O <sub>6</sub>
<b>Purity</b>	> 95% by hplc
<b>MW</b>	954.16 Da
<b>CAS Number</b>	673202-67-0
<b>Solubility (25°C)</b>	Soluble to 2 mg/ml in water
<b>Expiration Date</b>	12 months from date of receipt.