

Product Datasheet

TAT-Gap19 (orb1140646)

Description	Connexin43 hemichannel inhibitor; Cell penetrating peptides.	Image not found or type unknown
Form/Appearance	Freeze dried solid	Structure for TAT-Gap19
Storage	Store dry, frozen and in the dark	
Note	For research use only	
Protein Sequence	YGRKKRRQRRRKQIEIKKFK	
Formula	$C_{119}H_{212}N_{46}O_{26}$	
Purity	> 95% by hplc	
MW	2703.28 Da	
CAS Number	1507930-54-2	
Solubility (25°C)	Soluble to 1 mg/ml in water	
Expiration Date	12 months from date of receipt.	