

Product Datasheet

WaTx (orb1140661)

Description	Stabilizes the transient receptor potential cation channel 1 (TRPA1) in an active state; Peptides.
Form/Appearance	Freeze dried solid
Storage	Storage desiccated, frozen and in the dark
Note	For research use only
Protein Sequence	ASPQAKYCYEQCNVKNKVPFDQCYQMCSPLERS
Formula	$C_{164}H_{245}N_{45}O_{53}S_5$
Purity	> 95% by HPLC
MW	3855.3 Da
Solubility (25°C)	Soluble in DMSO as stock.
Expiration Date	12 months from date of receipt.