

Product Datasheet

Cusatuzumab biosimilar antibody (orb1140829)

Description	Cusatuzumab biosimilar antibody
Conjugation	Unconjugated
Preservatives	PBS, pH7.5
Concentration	1mg/ml
Storage	It is recommended to aliquote and store at -20°C for long-term storage. Avoid repeated freezing and thawing cycles.
Note	For research use only
Isotype	IgG1 Lambda
Purity	> 95% as determined by SDS-PAGE
Source	CHO Cells
CAS Number	1399546-01-0
Expiration Date	12 months from date of receipt.