

## Product Datasheet

### Belimumab biosimilar antibody (orb1141144)

<b>Description</b>	Belimumab biosimilar antibody
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, FC, Functional Studies, IF, IP, NeA
<b>Concentration</b>	1mg/ml
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1 Lambda
<b>Purity</b>	≥95%
<b>Expiration Date</b>	12 months from date of receipt.