

## Product Datasheet

# Rabbit anti-COVID-19/SARS-CoV-2 S1-NTD Polyclonal Antibody (orb1141354)

<b>Description</b>	Rabbit anti-COVID-19/SARS-CoV-2 S1-NTD Polyclonal Antibody
<b>Species/Host</b>	Rabbit
<b>Endotoxins</b>	< 0.1EU/ug determined by LAL method.
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Expiration Date</b>	12 months from date of receipt.