

Product Datasheet

Rabbit anti-COVID-19/SARS-CoV-2 S1-RBD Polyclonal Antibody (orb1141357)

Description	Rabbit anti-COVID-19/SARS-CoV-2 S1-RBD Polyclonal Antibody
Species/Host	Rabbit
Endotoxins	< 0.1EU/ug determined by LAL method.
Conjugation	Unconjugated
Tested Applications	ELISA
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.