

## Product Datasheet

### Recombinant Human ACE2 Protein (orb1141403)

<b>Description</b>	Recombinant Human ACE2 Protein
<b>Species/Host</b>	Mammal
<b>Note</b>	For research use only
<b>Purity</b>	≥ 90% as determined by SDS-PAGE
<b>Source</b>	Recombinant Human ACE2 was expressed from Mammalian cells expression system and the target gene encoding Ser19-Asp615 is expressed with C-His Tag
<b>Expiration Date</b>	6 months from date of receipt.