

Product Datasheet

Recombinant SARS-CoV-2/COVID-19 ORF8 Protein (orb1141412)

Description	Recombinant SARS-CoV-2/COVID-19 ORF8 Protein
Species/Host	E. coli
Note	For research use only
Purity	≥ 90% as determined by SDS-PAGE
Source	Recombinant SARS-CoV-2 ORF8 protein was expressed from E.coli expression system and the target gene encoding Phe16-Ile121 is expressed with N-His Tag
Expiration Date	6 months from date of receipt.