

Product Datasheet

Recombinant SARS-CoV-2/COVID-19 3C-like Proteinase (orb1141417)

Description	Recombinant SARS-CoV-2/COVID-19 3C-like Proteinase
Species/Host	E. coli
Note	For research use only
Purity	≥ 90% as determined by SDS-PAGE
Source	Recombinant SARS-CoV-2 3C-like Proteinase was expressed from E.coli expression system and the target gene encoding Ser1-Thr304 is expressed with N-GST Tag
Expiration Date	6 months from date of receipt.