

## Product Datasheet

# Recombinant SARS-CoV-2/COVID-19 3C-like Proteinase (orb1141419)

<b>Description</b>	Recombinant SARS-CoV-2/COVID-19 3C-like Proteinase
<b>Species/Host</b>	E. coli
<b>Note</b>	For research use only
<b>Purity</b>	≥ 90% as determined by SDS-PAGE
<b>Source</b>	Recombinant SARS-CoV-2 3C-like Proteinase was expressed from E.coli expression system and the target gene encoding Ser1-Gln306 is expressed with N-His Tag
<b>Expiration Date</b>	6 months from date of receipt.