

## Product Datasheet

# Recombinant SARS-CoV-2/COVID-19 NSP3 Protein (orb1141420)

<b>Description</b>	Recombinant SARS-CoV-2/COVID-19 NSP3 Protein
<b>Species/Host</b>	E. coli
<b>Note</b>	For research use only
<b>Purity</b>	≥ 90% as determined by SDS-PAGE
<b>Source</b>	Recombinant SARS-CoV-2 NSP3 was expressed from E.coli expression system and the target gene encoding Glu1024-Gln1198 is expressed with N-His Tag
<b>Expiration Date</b>	6 months from date of receipt.