

Product Datasheet

Human Interferon γ (IFN- γ) ELISA Kit (orb1141475)

Description	Human Interferon γ (IFN- γ) ELISA Kit
Reactivity	Human
Range	6.25 pg/mL-400 pg/mL
Tested Applications	ELISA
Target	IFNG
Note	For research use only
Sample Types	serum, cell culture supernates, urine, tissue homogenates, cerebrospinal fluid (CSF)
Uniprot ID	P01579
Sensitivity	1.56 pg/mL
Expiration Date	Please enquire.