

Product Datasheet

IRF7 (phospho-Ser471/472) antibody (Biotin) (orb114208)

Description	IRF7 (phospho-Ser471/472) antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	ELISA, ICC, IHC-P, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human IRF7 around the phosphorylation site of Ser471/472
Target	Phospho-IRF7 (Ser471 + Ser472)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q92985
Dilution Range	WB=1:500-2000,IHC-P=1:100-500,ICC=1:100,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.