

Product Datasheet

EGLN2 antibody (orb1144175)

Description	Mouse monoclonal antibody to EGLN2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human EGLN2(NP_542770) produced in E.coli.
Preservatives	PBS (PH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	43.5 kDa(Predicted)
Uniprot ID	A0A024R0R2
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.