

Product Datasheet

OXSR1 antibody (orb1144950)

Description	Mouse monoclonal antibody to OXSR1
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full-length protein expressed in 293T cell transfected with human OXSR1 expression vector
Preservatives	PBS (pH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	57.8 kDa(Predicted)
Uniprot ID	O95747
Dilution Range	WB 1:1000, IHC 1:150
Expiration Date	12 months from date of receipt.