



Product Datasheet

PCBP2/hnRNP E2 Antibody (monoclonal, 4B9C7) (orb1145856)


Description	PCBP2/hnRNP E2 Antibody (monoclonal, 4B9C7)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, WB
Immunogen	E.coli-derived human PCBP2/hnRNP E2 recombinant protein (Position: Q197-K276).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Isotype	Mouse IgG2a
Clonality	Monoclonal
MW	39 kDa
Uniprot ID	Q15366
Expiration Date	12 months from date of receipt.



IF analysis of PCBP2/hnRNP E2 using anti-PCBP2/hnRNP E2 antibody.



Western blot analysis of PCBP2/hnRNP E2 antibody.



Flow Cytometry analysis of PC-3 cells showing PCBP2/hnRNP E2 expression.