



## Product Datasheet

### EIF5 Antibody (monoclonal, 2F13C4) (orb1145868)


<b>Description</b>	EIF5 Antibody (monoclonal, 2F13C4)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC, ICC, IF, WB
<b>Immunogen</b>	E.coli-derived human EIF5 recombinant protein (Position: D11-D423).
<b>Form/Appearance</b>	Lyophilized
<b>Concentration</b>	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
<b>Isotype</b>	Mouse IgG2b
<b>Clonality</b>	Monoclonal
<b>MW</b>	49 kDa
<b>Uniprot ID</b>	<b>P55010</b>
<b>Expiration Date</b>	12 months from date of receipt.



IF analysis of EIF5 using anti-EIF5 anti...



Western blot analysis of EIF5 using anti...



Flow Cytometry analysis of SiHa cells us...