

Product Datasheet

Mouse pMAPT/pTAU(phosphorylated microtubule-associated protein tau) ELISA Kit (orb1146714)

Description The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

Reactivity Mouse

Range 15.63-1000 pg/mL

Note For research use only

Sample Types serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids

Sensitivity 9.38 pg/mL

Expiration Date Please enquire.

