

Product Datasheet

Guinea pig ADMA(Asymmetrical Dimethylarginine) ELISA Kit (orb1146781)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Guinea pig
Range	15.63-1000 ng/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Sensitivity	4.51 ng/mL
Expiration Date	Please enquire.

