

## Product Datasheet

# Rat TRPV1(Transient Receptor Potential Cation Channel Subfamily V, Member 1) ELISA Kit (orb1146838)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Rat
<b>Range</b>	0.16-10 ng/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	Tissue homogenates, cell lysates and other biological fluids.
<b>Uniprot ID</b>	<b>035433</b>
<b>Sensitivity</b>	0.059 ng/mL
<b>Expiration Date</b>	Please enquire.

