

Product Datasheet

Rat SLC(Secondary Lymphoid Tissue Chemokine) ELISA Kit (orb1146840)

| | |
|------------------------|---|
| Description | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| Reactivity | Rat |
| Range | 7.82-500 pg/mL |
| Note | For research use only |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids |
| Sensitivity | 3.2 pg/mL |
| Expiration Date | Please enquire. |

