

## Product Datasheet

# Human TrkA/NTRK1(High affinity nerve growth factor receptor) ELISA Kit (orb1146845)

|                        |   |
|------------------------|---|
| <b>Description</b>     | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| <b>Reactivity</b>      | Human   |
| <b>Range</b>           | 78.13-5000 pg/mL  |
| <b>Note</b>            | For research use only   |
| <b>Sample Types</b>    | serum, plasma and other biological fluids   |
| <b>Uniprot ID</b>      | <b>P04629</b>   |
| <b>Sensitivity</b>     | 31 pg/mL  |
| <b>Expiration Date</b> | Please enquire.   |

