

Product Datasheet

Human SCUBE1(Signal peptide, CUB and EGF-like domain-containing protein 1) ELISA Kit (orb1146867)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	0.16-10 ng/mL
Note	For research use only
Sample Types	Tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q8IWY4
Sensitivity	0.094 ng/mL
Expiration Date	Please enquire.

