

Product Datasheet

Human BNIP3L(BCL2/adenovirus E1B 19 kDa protein-interacting protein 3-like) ELISA Kit (orb1146871)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	0.16-10 ng/mL
Note	For research use only
Sample Types	Tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	O60238
Sensitivity	0.094 ng/mL
Expiration Date	Please enquire.

