

## Product Datasheet

### Sheep HP(Haptoglobin) ELISA Kit (orb1146881)

|                        |   |
|------------------------|---|
| <b>Description</b>     | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| <b>Reactivity</b>      | Sheep   |
| <b>Range</b>           | 15.63-1000 ng/mL  |
| <b>Note</b>            | For research use only   |
| <b>Sample Types</b>    | serum, plasma and other biological fluids   |
| <b>Sensitivity</b>     | 5.9 ng/mL   |
| <b>Expiration Date</b> | Please enquire.   |

