

## Product Datasheet

# 10-HETE(10-Hydroxyeicosatetraenoic Acid) ELISA Kit (orb1146884)

<b>Description</b>	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
<b>Reactivity</b>	All
<b>Range</b>	12.5-800 ng/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma and other biological fluids
<b>Sensitivity</b>	3.95 ng/mL
<b>Expiration Date</b>	Please enquire.

