

Product Datasheet

Human sTNFRSF1A(Soluble Tumor Necrosis Factor Receptor Superfamily, Member 1A) ELISA Kit (orb1146887)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	15.63-1000 pg/mL
Note	For research use only
Sample Types	serum, plasma, urine, cell lysates and other biological fluids
Sensitivity	6.5 pg/mL
Expiration Date	Please enquire.

