

Product Datasheet

Pig GnRH(Gonadotropin Releasing Hormone) ELISA Kit (orb1146939)

| | |
|------------------------|---|
| Description | This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro... |
| Reactivity | Porcine |
| Range | 15.63-1000 pg/mL |
| Note | For research use only |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids |
| Uniprot ID | P49921 |
| Sensitivity | 5.25 pg/mL |
| Expiration Date | Please enquire. |

