

## Product Datasheet

# Rat GSTs(Glutathione S Transferases) ELISA Kit (orb1146940)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Rat
<b>Range</b>	0.16-10 ng/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
<b>Uniprot ID</b>	<a href="#">O35543</a>
<b>Sensitivity</b>	0.056 ng/mL
<b>Expiration Date</b>	Please enquire.

