

## Product Datasheet

# Rat $\tau$ (phospho Tau Protein) ELISA Kit (orb1146942)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Rat
<b>Range</b>	15.63-1000 pg/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
<b>Sensitivity</b>	5.7 pg/mL
<b>Expiration Date</b>	Please enquire.

