

## Product Datasheet

# Mouse cADPRH(Cyclic ADP Ribose Hydrolase) ELISA Kit (orb1146953)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Mouse
<b>Range</b>	1.57-100 ng/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	tissue homogenates, cell lysates and other biological fluids
<b>Uniprot ID</b>	<b>P56528</b>
<b>Sensitivity</b>	0.55 ng/mL
<b>Expiration Date</b>	Please enquire.

