
Product Datasheet

Human PDGF R beta / CD140b Protein (orb1146967)

Description

Human PDGF R beta, Fc Tag is

(HEK293). It contains AA Leu 33 - Phe 53...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

PDGF R beta

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 83.1 kDa. The protein migrates as 100-120 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

NP_002600.1

Purity

95%

Source

Human PDGF R beta, Fc Tag (orb1146967) is expressed from human 293 cells (HEK293). It contains AA Leu 33 - Phe 530 (Accession # P09619-1).

MW

83.1 kDa

Uniprot ID[P09619](#)**NCBI**[NP_002600.1](#)**Expiration Date**

6 months from date of receipt.