

Product Datasheet

Human IL-5 R alpha / CD125 Protein (orb1146968)

Description

Human IL-5 R alpha, Fc Tag is expressed from

Reactivity

Trp 3...

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

IL-5 R alpha

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

N-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the N-terminus. The protein has a calculated MW of 63.2 kDa. The protein migrates as 70-90 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

95%

Source


Human IL-5 R alpha, Fc Tag (orb1146968) is expressed from human 293 cells (HEK293). It contains AA Asp 21 - Trp 342 (Accession # Q01344-1).

MW

63.2 kDa

Uniprot ID
[Q01344](#)
NCBI
[NP_000555.2](#)
Expiration Date

6 months from date of receipt.



Immobilized
Human IL-5,
His Tag
(Cat. No...