

Product Datasheet

Cynomolgus IL-1 RAcP / IL-1 R3 Protein (orb1146971)

Description

Cynomolgus IL-1 RAcP, His Tag is expressed from

Reactivity	Gl...
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IL-1 RAcP
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 41.0 kDa. The protein migrates as 55-65 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	Cynomolgus IL-1 RAcP, His Tag (orb1146971) is expressed from human 293 cells (HEK293). It contains AA Ser 21 - Glu 359 (Accession # G7NYP7-1).
MW	41.0 kDa
Uniprot ID	G7NYP7
Expiration Date	6 months from date of receipt.

 Cynomolgus IL-1 RAcP His Tag EISA
 1.1 µg of human Anti-IL1RAP human IgG1 per well

 Immobilized
 Human
 Anti-IL1RAP
 Human
 IgG1...