


Product Datasheet

Human IGFBP-7 Protein (orb1146981)

Description	Human IGFBP-7, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Asp 30 - Leu 282 (A...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IGFBP-7
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 52.6 kDa. The protein migrates as 55-66 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	90%
Source	Human IGFBP-7, Fc Tag (orb1146981) is expressed from human 293 cells (HEK293). It contains AA Asp 30 - Leu 282 (Accession # Q16270-1).
MW	52.6 kDa
Uniprot ID	Q16270
NCBI	NP_001544.1
Expiration Date	6 months from date of receipt.



Human
IGFBP-7
Protein, Fc
Tag (Cat.
No. ...