


Product Datasheet

Human IL-20 R alpha Protein (orb1146994)

Description	Human IL-20 R alpha, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Val 30 - Lys ...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IL-20 R alpha
Preservatives	PBS, pH 7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 51.7 kDa. The protein migrates as 65-70 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	Human IL-20 R alpha, Fc Tag (orb1146994) is expressed from human 293 cells (HEK293). It contains AA Val 30 - Lys 250 (Accession # Q9UHF4-1).
MW	51.7 kDa
Uniprot ID	Q9UHF4
NCBI	NP_055247.3
Expiration Date	6 months from date of receipt.



Human IL-20, His Tag immobilized on CM5 ...