

Product Datasheet

Rat DLL3 Protein (orb1147001)

Description	Rat DLL3, His Tag is expressed from human 293 cells (HEK293). It contains AA Ala 25 - Ala 488 (Acces...
Reactivity	Rat
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	DLL3
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 50.1 kDa. The protein migrates as 50-60 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_446118.1
Purity	90%
Source	Rat DLL3, His Tag (orb1147001) is expressed from human 293 cells (HEK293). It contains AA Ala 25 - Ala 488 (Accession # O88671-1).
MW	50.1 kDa
Uniprot ID	O88671
NCBI	NP_446118.1
Expiration Date	6 months from date of receipt.