

Product Datasheet

Human BTLA (31-150) Protein (orb1147004)

Description	Human BTLA (31-150), His Tag is expressed from human 293 cells (HEK293). It contains AA Lys 31 - Ser...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	BTLA
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 15.7 kDa. The protein migrates as 27-34 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	AAP44003.1
Purity	90%
Source	Human BTLA (31-150), His Tag (orb1147004) is expressed from human 293 cells (HEK293). It contains AA Lys 31 - Ser 150 (Accession # AAP44003.1).
MW	15.7 kDa
NCBI	AAP44003.1
Expiration Date	6 months from date of receipt.